June 2024

QIAGEN Digital Insights Admin Tool Quick Start Guide



Contents

Introduction
Access to QDIAT
How to use QDIAT4
Administrators
Add New Users to Licensing Group5
Users in License Group8
Remove Users from License Group to End Access8
Sessions10
Viewing who is currently logged into your license group10
Terminating an active session10
Viewing failed sessions10
Usage Monitoring11
License Usage12
User Usage for Sample Upload12
SSO Configuration13
Prerequisites13
SSO Config: step 114
SSO Config: step 215
SSO Config: step 316
How to navigate between the SSO configuration steps
Changes to current support process for authentication related issues17

Introduction

QIAGEN Digital Insights Admin Tool (QDIAT) is a platform that allows you to:

- 1. Manage your users and monitor your license(s)/sessions/usage
- 2. Single Sign-on (SSO) configuration

This platform manages the access to the following products

- HGMD Human Genetic Mutation Database -<u>https://my.qiagendigitalinsights.com/bbp/</u>
- HSMD Human Somatic Mutation Database -<u>https://my.qiagendigitalinsights.com/bbp/</u>
- COSMIC Catalogue of Somatic Mutations in Cancer -<u>https://my.qiagendigitalinsights.com/bbp/</u>
- PGMD PharmacoGenomic Mutation Database -<u>https://my.qiagendigitalinsights.com/bbp/</u>
- ANNOVAR https://my.qiagendigitalinsights.com/bbp/
- Genome Trax[™] <u>https://my.qiagendigitalinsights.com/bbp/</u>
- QCI QIAGEN Clinical Insights https://variants.ingenuity.com/qci
- IPA Ingenuity Pathway Analysis https://analysis.ingenuity.com/pa

QDIAT allows you to look up which licenses users have within your organization, control access to products by adding or removing users, look up who is currently logged in, and monitor when licenses expire.

Access to QDIAT

If you are a license coordinator for any of these products, log in at <u>https://apps.ingenuity.com/iat</u> with your username and password.

If you are the license coordinator and cannot access QDIAT please contact customer support via <u>TS-Bioinformatics@qiagen.com</u> and request a QDIAT license. Please note, a license to your

product can be managed by more than one administrator. If you wish to add/remove administrators, you must contact Customer Support. Administrators do not need to be members of the License Group (individual users that have access to a licensed product).

How to use QDIAT

License Administrator/License Coordinator: this is the administrator at the institution who is authorized to make changes to the licenses/users.

Once logged in, QDIAT displays on the home page the list of active licensing groups. From here you can administer more than one product license. To do this, you will need to choose which license you want to access after you log into QDIAT. In the example shown below, the license administrator can manage two different licenses, and must select one before continuing.

QIAGEN Digital Insights Admin Tool	Home	SSO configuration	mihaela.suteu@qiagen.com	Logout
Home				
My Licensing Groups				
QIAGEN - IPA - Internal				
QIAGEN - QCI - Internal				
QIAGEN - MyQDI - Internal				

Once a license group is selected, the top menu will display the following tabs: Users and "QIAGEN Digital Insights" with the following submenu: Licenses, Sessions, Usage, User Groups and User Features.

QIAGEN Digital Insights Admin Tool Users	QIAGEN Clinical Insight Interpret	mihaela.suteu@qiagen.com	Logout
	Licenses Sessions Usage User Groups User Features		
Home / QIAGEN - QCI - Internal / QIAGEN Clinical I	sight Interpret / Licenses		
Licenses			

User Management

In the Users tab for QDIAT you will find 3 sections; (1) Administrators, (2) Add New Users to Licensing Group, and (3) Users in Licensing Group. These sections are described in more detail below.

	N Digital Insights Admin Tool	Users	About			mihaela.suteu	@qiagen.com
-lome	QIAGEN - QCI - Internal		All users listed in the Users tab belong to the license groups and have access to all applications shown.				
Add U	Isers to Licensing Group (2)				Administrators 1		
Sear	ch Users Create Account				Username	Name	Admin Type
Search	and add existing accounts to the licer	nsing group.			christine.bueno@qiagen.com	Christine Bueno	SSO
User	account email address or addresses				cristian.balint@qiagen.com	Cristian Balint	SSO
Enter or	ne or more email addresses or search	on one parti	al address.		mihaela.suteu@qiagen.com	Mihaela Suteu	License
Searc	ch				roslen.anacleto@qiagen.com	Roslen Anacleto	SSO
					sachin.pattanshetty@qiagen.com	Sachin Pattanshetty	SSO
					stuart.tugendreich@qiagen.com	Stuart Tugendreich	SSO
show o	in Licensing Group 3 only users with SSO Export users 10 entries	as: XLS XL	SX		stuart.tugendreich@qiagen.com	Stuart Tugendreich	SSO Search Date
how o	only users with SSO Export users	as: XLS XL	SX		Search Name		
how o	only users with SSO Export users	as: XLS XL		Nar	Search Name		Search Date
how o	nly users with SSO Export users 10 entries Search Username Username	as: XLS XL		Nar	Search Name		Search Date
Show o	nhy users with SSO Export users 10 entries Search Username Username aaron.mandel@qlagen.com	as: XLS XL		Aar Aar	Search Name me on Mandel		Search Date Added Date May 19 2020




Administrators

Displays a list of administrators that can manage the access of users to the product. This list is managed by Customer Support. If you wish to add/remove administrators, you must contact Customer Support via <u>IS-Bioinformatics@giagen.com</u>.



Add New Users to Licensing Group

It is possible that users within your organization possess licenses for other QIAGEN Digital Insights products and are already registered within the QIAGEN licensing system. If this is the case, you can simply search in the "Add User to Group" section to see if the user has an account, as shown below. Note that in the screenshot below the search is performed using the organizational email address.



If the user exists within QIAGEN's licensing system a window will pop-up displaying the user details. You can then select the user to be added to this license group. Click on "Save Changes" to finalize the process and user will receive a "Welcome" email with a link to access this product.

Add Users To Li	censing Group	
Select the users the	at will be added to the licensing group QIAGEN Internal - MyQDI	
Show 10 🗸	untries	
S.	Username	Name 🔶
	dominic.john@qiagen.com	Dominic John
Showing 1 to 1 of	1 entries	Previous 1 Next
		Cancel Save Changes

If the user is not found, you can create a user account for that user. To add a user to your license group, click on "Create Account", fill in the required information, and click "Add to License".

Add Users to Licensing Group	
Search Users Create Account	
Create new account and add it to the licensing group.	
Institutional Email Address	
Enter Institutional Email Address	
First Name	
Enter First Name	
Last Name	
Enter Last Name	
Institution or Company Name	
The license does not have a domain or the user's domain is not configured for SSO	h
Enable SSO	
Add to License	Reset

When entering the Institution or Company name, suggestions will appear via autocomplete. Please try to use one of these names if applicable. If the institution or company is not listed, then type in the name.

Checking the *Enable SSO* checkbox will enable the account to login via the institution's Identity Provider. If the checkbox is in read-only mode, this means there is no SSO configured for your domain's institution.

After a moment, you will see a confirmation that the user was created and added to your license group. The user will then receive a "Welcome" email to activate account, set password and link to access this product.

Add Users to Lice	ensing Group
Search Users	Create Account
The new accoun group.	t was successfully created and added to the licensing $\qquad \times$
Create new accoun	nt and add it to the licensing group. I Address
Enter Institutional	I Email Address
First Name	
Enter First Name	
L (N	
Last Name	



Users in License Group

This section displays users in your license group for your product. The table is searchable by email, name, and the date they were added to the license. A filter can be applied to view only the users in the license group having SSO enabled.

If an account has SSO enabled, a key icon will be displayed near the username.

The list of users can be exported in a xls or xlsx file format.

Users	in Licensing Group 3		
Show	only users with SSO Export users as: XLS XLSX		
Show	10 v entries		
	Search Usemame	Search Name	Search Date
•	Username	Name ÷	Added Date
0	▶ aarthi.goverdhan@qiagen.com	Aarthi Goverdhan	Nov 04 2021
0		Aavishkar Ghosh	Sep 20 2022
0	ada.cean@contractor.qiagen.com	Ada Cean	Nov 19 2020



Remove Users from License Group to End Access.

To remove users from your license, first check the box to the left of the member(s) you wish to remove and then click on the "Remove From Licensing Group" button.

Users in	Licensing Group		
Show	10 v entries		Remove From Licensing Group
	Beate	Search Name	Search Date
ð	Username	Name ¢	Added Date
Ø	beate.litzenburger@qiagen.com	Beate Litzenburger	Aug 13 2020
Showin	115 of 1 entries (filtered from 129 total entries)		Previous 1 Next

Once you remove a user from a specific licensed product, any active session for that person is immediately terminated and he or she will not be able to access this specific product.

License Management

The Licenses tab displays the start and end date of active license(s) for the license group.

IAGEN Digital Insight	s Admin Tool Users QIAGEN	Olinical Insight Interpret					mihaela.suteu@qiagen.com	Logo
	Licenses	Sessions Usage User	Groups					
Home / QIAGEN - QCI - I	Internal / QIAGEN® Clinical Insight Inte	rpret / Licenses						
Licenses								
Show 10 v entries								
License ID	Application Name	Activation Date	End Date	IP addresses 0	Domains 😧 🔶	2FA Duration 🔮 🔶	MPC Duration	4
5908912	QIAGEN Clinical Insight Interpret	Sep 01 2020	Aug 29 2030	-	-	60 days	365 days	
5908913	QIAGEN Clinical Insight Interpret	Sep 01 2020	Aug 29 2030	-	-	60 days	365 days	
Showing 1 to 2 of 2 entries	5						Previous 1	Next

If the license restricts usage from specific IP addresses, those will be listed. Then all sessions on this license must originate from these IP addresses. This is the institutional/building IP address as seen by the rest of the world. This is not the IP address on the user's computer.

Please note, only users from the specified domains can use your license, which are the email domains belonging to the institution that bought the license (e.g. qiagen.com).

Security settings defined on the license are also displayed here. These can be custom defined values or default ones provided by the system.

The two factor authentication duration column lists the duration for which a 2FA token is valid for, it applies to users that are not on SSO. If the license group contains both users that are on SSO login and users that are not on SSO, then the 2FA setting will be displayed as a combination.

The mandatory password change duration column lists the duration for which a password is valid for, after which the system will request to be changed with a new one. It applies to users that are not on SSO. If the license group contains both users that are on SSO login and users that are not on SSO, then the MPC setting will be displayed as a combination.

Sessions

Viewing who is currently logged into your license group.

Go to the *Live Sessions* section, which provides a list of users who are currently logged into your license. You can see how long they have been using the license in the current session and you can terminate that session if necessary (see below). Note that all times listed are in the PST time zone.

GEN Digital Insights Admin Too	Users C	IAGEN Clinica	I Insight Interpr	er			minaela.su	teu@qiagen.com	Lo
	L	icenses Se	ssions Usag	е					
ome / QIAGEN - QCI - Internal / QIAG	EN Clinical Insigh	nt Interpret / Se	essions						
cense Details									
cense ID:	956531								
lax. Concurrent Users:	2147483	647							
how 10 v entries									
Jsername	Session Id	Start Time $^{\diamondsuit}$	Last Request	Duration	Inactivity Timeout	IP Address	Geo Location	Action	
mihaela.suteu@qiagen.com	18710129	Aug 31 2023, 16:22:02	Aug 31 2023, 18:17:56	1h 55m	1h 0m	10.32.45.46	Rotterdam,South Holland,Netherlands	Terminate	
onel.marincas@qiagen.com	18710129	Aug 31 2023, 16:22:02	Aug 31 2023, 18:17:56	1h 55m	1h 0m	10.32.45.46	Rotterdam,South Holland,Netherlands	Terminate	

Terminating an active session

In the Live Sessions section, you can click the terminate button to the right of the user whose session you wish to terminate. Please use this option only in urgent situations, as it will potentially interrupt the user's work. We recommend that you contact the user before terminating their session.

Viewing failed sessions

The Failed Sessions section provides a list of users who could not use the application. It states the cause why the user could not access the application.

ailed Sessions						
now 10 v entries						
Username	∲ User Id	Global Session Id	Application [‡]	Timestamp 🔻	Failure Cause	IP Address
tasim.banu@qiagen.com	432251	515373287	IAT	08/09/2021 04:07:03	NO_VALID_LICENSE_FOR_APP	106.200.247.139
haijin.lei@qiagen.com	340256	514426663	IAT	08/05/2021 21:51:00	NO_VALID_LICENSE_FOR_APP	58.63.130.168
haijin.lei@qiagen.com	340256	514162058	IAT	08/05/2021 05:37:40	NO_VALID_LICENSE_FOR_APP	58.33.58.146
haijin.lei@qiagen.com	340256	0	INGSSO	08/05/2021 05:37:29	PASSWORD_INVALID	58.33.58.146

Usage Monitoring

The usage tab allows users to review the sample upload limit and sample consumption associated with your license. This tab Is divided into 2 sections that will be explained in more detail below.

QIAGEN Clinical Insight Interpret Aug 31 2020 Aug 28 2030		Show 10 v en
eature Name	Feature Description	Feature Current Usage / Limit
CIT Upload	QCIT Upload	19 / 10000
Jpload (2)		Previous 1 Ne
<u> </u>		Show 10 v ent
sername 🔶		Individual Usage
an-noel.billaud@qiagen.com	Jean-Noel Billaud	Individual Usage 0 / 10000
an-noel billaud@qiagen.com aron.mandel.qcit@qiagen.com	Jean-Noel Billaud aaron mandel_qcit	0 / 10000 0 / 10000
an-noel. billaud@qiagen.com aron.mandel.qcit@qiagen.com s_selenium_qcit_autouser@qiagen.com	Jean-Noel Billaud aaron mandei_qcit VCS Selenium QCI-T Auto User	0 / 10000 0 / 10000 0 / 10000
an-noel. billaud@qiagen.com aron.mandel.qcit@qiagen.com :s_selenium_qcit_autouser@qiagen.com m@demo.com	Jean-Noel Billaud aaron mandel_qcit VCS Selenium QCI-T Auto User Tim Demo	Individual Usage 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000
an-noel.billaud@qiagen.com aron.mandel.qcit@qiagen.com :s_selenium_qcit_autouser@qiagen.com m@demo.com ong@demo.com	Jean-Noel Billaud aaron mandel_qcit VCS Selenium QCI-T Auto User Tim Demo Dong Demo	Individual Usage 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000
an-noel.billaud@qiagen.com aron mandel.qcit[@qiagen.com :s_selenium_qcit_autouser@qiagen.com m@demo.com ong@demo.com mes@demo.com	Jean-Noel Billaud aaron mandel_qcit VCS Selenium QCI-T Auto User Tim Demo	Individual Usage 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000
an-noel billaud@qiagen.com aron mandel.qcit@qiagen.com is_selenium_qcit_autouser@qiagen.com m@demo.com ong@demo.com mes@demo.com nawn@demo.com	Jean-Noel Billaud aaron mandel_qcit VCS Selenium QCI-T Auto User Tim Demo Dong Demo James Demo Shawn Demo	Individual Usage 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000
an-noel.billaud@qiagen.com aron mandel.qcit[@qiagen.com :s_selenium_qcit_autouser@qiagen.com m@demo.com ong@demo.com mes@demo.com	Jean-Noel Billaud aaron mandel_qcit VCS Selenium QCI-T Auto User Tim Demo Dong Demo James Demo	Individual Usage 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000 0 / 10000
	eature Name ¢ CIT Upload	eature Name Feature Description CIT Upload CIT Upload



License Usage

This section shows you the activation and end date of your QCI Interpret Translational license. Your QCI Interpret Translational license consumption limit and usage is noted in the table above. In the example, the Group has a limit of 10000 samples of which 19 samples have been consumed by the Group.



User Usage for Sample Upload

This table displays the number of samples consumed per user within the Group. From this table you can see the individual consumption rate.

SSO Configuration

The SSO Configuration tab allows the SSO Administrator to configure SSO for his/her company and to manage users by enabling/disabling the SSO process for them.

QIAGEN Digital Insights Admin Tool	Home	SSO configuration	mihaela.suteu@qiagen.com	Logout
		W		
Home				
My Licensing Groups				
QIAGEN - IPA - Internal				
QIAGEN - QCI - Internal				
QIAGEN - MyQDI - Internal				

SSO Administrator: this is typically the buyer's IT person who configures and manages the institutions single sign-on with QDI's product offerings. To configure SSO for your institution reach out to <u>TS-Bioinformatics@qiagen.com</u>. We support multiple SSO Administrators.

Qiagen Digital Insights is an identity provider that offers authentication with security measures such as complex passwords, password expiration, 2 factor authentication. One user account in our identity store can access all QDI products as listed above in the Introduction.

Customers have the option to use their own identity provider for authentication through QDI's Single Sign-on (SSO). SSO allows you to login using your company credentials. Qiagen's SSO is based on SAML 2.0, more specifically, the Service Provider-Initiated SAML Flow.

Once SSO is configured, the SSO Administrator is responsible for all communication and authentication related issues/queries for their users.

Prerequisites

- Send a request to <u>TS-Bioinformatics@giagen.com</u>
 - \circ ~ with first name, last name, and email ID of the SSO Administrator
 - list of domains for SSO

An SSO administrator will be created by the QDI customer support team and in order to gain access to the SSO configuration tab, a special license will be assigned.

Creating a new configuration

On the right side panel you can find detailed instructions on how to create and setup a new SSO configuration for your company.

On the left side panel, in the My SSO Configs section you will find a list of your current SSO configurations. You can click on the Add new SSO config to create a new one.


SSO configuration	
My SSO Configs	SSO instructions
lon.com	Overview
keycloak.com	Single sign-on allows you to login using your company credentials. QIAGEN single sign-on (SSO) is based on SAML 2.0. QIAGEN can map attributes to provision a user to different group with feature controls.
glagen.com contractor.glagen.com	CALCEN acts as the Service For Vorder (SP) and It's capable of offering automatic user provisioning. In case auto-provisioning is enabled, then you do not need to register as a user in QIAGEN. Once QIAGEN receives a SAML response from the Identity Provider (IdP), QIAGEN tocks if this user exists. If the user does not exist, QIAGEN creates a user account automatically with the received name ID.
Administrators	This article covers: Configuring SSO Enabling SSO
Username 🔺 Name 💠 Admin Type	You can download QIAGEN service provider public certificate used for SAML signing from here
diana.stanila@qiagen.com Diana Stanila SSO	Configuring SSO
gabriel.padurean@qiagen.com Gabriel Padurean SSO	First, enter your SSO information
ionel.marincas@qiagen.com ionel Marincas License, SSO	Domain name(s): <domains do="" for="" sso="" to="" want="" which="" you=""> Identity provider id: <id entitydescriptor="" of=""></id></domains>
suteum@qiagen.com Miha Suteu SSO	 Identity provider metadata xml: <identity metadata="" provider="" xml=""> Note: The metadata xml should be provided to you from</identity>
	the IT department of your company Email attribute name: <name attribute="" email="" of=""> Note: The name of attribute should be provided to you from the IT</name>
	department of your company First attribute name: <name attribute="" first="" name="" of="">. Note: The name of attribute should be provided to you from the IT</name>
	 First attribute name: Note: The name of attribute should be provided to you from the T department of your company
	 Last attribute name: <name attribute="" last="" name="" of=""> Note: The name of attribute should be provided to you from the IT department of your company</name>
	Enabling SSO
	Second, download our metadata and enable SSO
	Note Once configured you'll have to download the SP metadata XML. Enable SSO
	Enabling SSO users
	Third, enable SSO login for existing users
	Note By default existing users will not be able to login with SSO. They have to be changed

There is also a table that lists all the administrators defined on the group.

There are three steps for setting up a new SSO configuration.

SSO Config: step 1

The configuration step requires you to choose one or multiple domain(s). The available domains are the ones setup by customer support at the time the SSO administrator was created. If you cannot find the domain you are looking for please reach out to customer support for any changes. Please observe the instructions on the right side panel for the information required in the form.



If the Automatic user provisioning is checked this means that all new users signing in to Qiagen via your company's identity provider will gain access to the Qiagen products by having accounts automatically created in the Qiagen system and added to your license group.

SSO Config: step 2

The 2nd step allows you to preview the configuration, choose the license group for auto provisioning and enable the SSO configuration. After selecting the domain(s) a service provider metadata URL will get displayed and it will allow for the QIAGEN metadata data file to be downloaded locally. This is also the step from which an existing SSO configuration can be disabled by clicking on the *Disable* button. A right side panel is also available with more information on this.

QIAGEN Digital Insights Admin Tool Home SSO configuration	Users	ioneLmarincas@qiagen.com Logou					
SSO configuration / demodomain.com							
(1) Configure SSO	Edit SS	3 Do configuration Enable SSO users					
View SSO configuration for domain(s)		SSO instructions					
Domain Name(s)							
demodomain.com		Overview					
Download QIAGEN service provider metadata		Single sign-on allows you to login using your company credentials. QIAGEN single sign-on (SSO) is based on SAML 2.0. QIAGEN can map attributes to provision a user to different group with feature controls.					
Entity Id		QIAGEN acts as the Service Provider (SP), and it's capable of offering automatic user provisioning. In case auto-provisioning is enabled, then you do not need to register as a user in QIAGEN. Once QIAGEN receives a SAML response from the Identity Provider					
demodomain.com		(dP), QIAGEN checks if this user exists. If the user does not exist, QIAGEN receives a seven response iron use network (dP), QIAGEN checks if this user exists. If the user does not exist, QIAGEN creates a user account automatically with the received name ID.					
xml version="1.0" encoding="UTF-8" standalone="yes"?		This article covers:					
<entitiesdescriptor< td=""><td></td><td>Configuring SSO</td></entitiesdescriptor<>		Configuring SSO					
xmlns="urn:oasis:names:tc:SAML:2.0:metadata" xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"		Enabling SSO You can download QIAGEN service provider public certificate used for SAML signing from here					
Name="http://ingenuity.com" veisebemal.eestian="umoosis:namesta:SAML-2.0;metadate.este.ea	mi cohomo motodato 3 D vod //	-					
Choose file No file chosen		Configuring SSO					
Email attribute name		First, enter your SSO information					
email		Domain name(s): <domains do="" for="" sso="" to="" want="" which="" you=""> Identity provider id: <id entitydescriptor="" of=""></id></domains>					
Enter name of attribute for email address sent in SAML response.		 Identity provider metadata xml: Identity provider metadata xml> Note: The metadata xml should be provided to you from the IT department of your company 					
First name attribute name		Email attribute name: <name attribute="" email="" of=""> Note: The name of attribute should be provided to you from the IT</name>					
firstname		department of your company First attribute name: <name attribute="" first="" name="" of=""> Note: The name of attribute should be provided to you from the IT</name> 					
Enter name of attribute for first name sent in SAML response.		department of your company Last attribute name: <name attribute="" last="" name="" of=""> Note: The name of attribute should be provided to you from the IT</name> 					
Last name attribute name		department of your company					
lastname		Enabling SSO					
Enter name of attribute for last name sent in SAML response.		Second, download our metadata and enable SSO					
Automatic user provisioning		Note Once configured you'll have to download the SP metadata XML.					
When this will be checked users will be automatically assigned to unlimited si	te licenses.	Enable SSO					
Show 10 V entries		Enabling SSO users					
License group name	Application name	Third, enable SSO login for existing users Note By default existing users will not be able to login with SSO. They have to be changed 					
	agenuity Views and Reports						
Demo Findings Internal Group In	ngenuity Views and Reports						
SSO Config for Demo	ngenuity Views and Reports						
Showing 1 to 3 of 3 entries	Previous 1 Next						
Save							
Save							

SSO Config: step 3

The 3rd step allows you to move users between two lists: one that is SSO enabled – meaning that users will sign in using your company identity provider and the other which is SSO disabled meaning that users will continue to sign in via Qiagen system.

In order to do this, you need to click on the checkbox next to the username and then use the arrows highlighted below to submit the action. The left and right side lists will refresh automatically after this.

IAGEN Digital Insights Admin Tool	Home SSO	configuration Users							demossoadmin@dem	odomain.com	Logo
SSO configuration / demodomain.com											
Configur			2 Disable SSO configuration	n		\leftarrow		Enab	3 le SSO users	•	
Users with demodomain.com in our s	ystem								•••		
Users with demodomain.com available fo	r SSO	Search:			I demodomain.co	m having SSO	enabled		Search:		
Show 10 V entries		Search:		Show 10	✓ entries				Search:		
Username	Full Name	Active licenses	÷	0	Username		Full Name	÷	Active licenses		÷
demossoadmin@demodomain.com	Demo SSOAdmin	QDI Admin Tool, User Manager Application	>				No data availa	ble in table			
foe.black@demodomain.com	Joe Black	Ingenuity Pathway Analysis, Ingenuity Views and Reports	<	Showing 0	to 0 of 0 entries					Previous	Next
mary.white@demodomain.com	Mary White	Ingenuity Pathway Analysis, Ingenuity Views and Reports									
to t	Tom Jones	Ingenuity Pathway Analysis, Ingenuity Views and Reports									
Showing 1 4 oventries		Previous 1 Next	t								

How to navigate between the SSO configuration steps

Use the arrows between the steps to move from one step to the other.

QIAGEN Digital Insights Admin Tool	SSO configuration	Users		demossoadmin@demodomain.com	Logout
SSO configuration / demodomain.com					
(1) Configure SSO		2 Disable SSO configuration	\leftarrow	3 Enable SSO users	
Users with demodomain.com in our sys	stem		W		

Changes to current support process for authentication related issues

When SSO is enabled, the SSO admin is responsible for all authentication related communication to their users.

Once QIAGEN assigns SSO admin for their domains/institution, QIAGEN will no longer reply/support authentication questions from those domains/institutions.